

Series 130 - Bioprocess Equipment Nozzle

Series BN Leistung BPE nozzle has a stainless steel body with tank interface surface, a gasket, a hygienic adapter and screwed fitting for adapter attachment. BPE nozzle is designed to be used as an aseptic connection to tank or vessel. Nozzle can be directly welded to the tank at the tank interface surface. Hygienic adapter is provided to connect with the flow-line. It is available in various end connections sizes and can be butt-welded directly to the tube. The gasket seal is provided to ensure completely safe and hygienic sealing of nozzle fittings.

- Excellent aseptic design
- Fully machined from SS 316L
- Unique screwed fittings
- Various internal surface finish available



Order Sample:

LE	130	A	20	E	T	9
1	2	3	4	5	6	7

Order Codes:

1	Make	Code
	Leistung	LE
2	Series	Code
	Bioprocess Equipment Nozzle	130
3	Connection std	Code
	ASME BPE	A
	BS 4825	B
	EN ISO 1127/ISO 4200	E
4	Connection Size	Code
	1/2"	15
	3/4"	20
	1"	25
	1 1/2"	40
	2"	50
	2 1/2"	65
	3"	80
	4"	100
	5"	125
	6"	150

5	Surface Finish	Code
	Mechanical	M
	Electropolish	E
6	Type of End Connection	Code
	Glass Fitment	G
	V - Band (Hygienic clamp)	T
7	Certificate	Code
	ISO 9001	9